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We are glad that you have chosen to join the thousands of satisfied users of WHITE ANGEL professional teeth whitening products.

Instructions from Your Dentist

Receive personal trays and consult your dentist regarding the use of the home whitening kit.

1. Brushing Teeth

Brush your teeth without toothpaste before starting the whitening process.

2. Applying the Whitening Gel

Follow your dentist's instructions on how to place the gel inside the personal whitening trays. Use no more than 1/3 of the syringe.

3. Inserting the Trays on Your Teeth

Insert the upper tray on the upper teeth and the lower tray on the lower teeth. Use a toothbrush to clean any excess gel that may spill onto the gums.

4. Waiting Period

For home whitening using WHITE ANGEL HOME KIT with a 20% , leave the trays in place for 2-4 hours per day. Continue the treatment for two to three weeks.

5. Finishing

After the time has passed, gently remove the trays and clean your mouth with a toothbrush and water. Then clean the trays to keep them ready for the next use.

Additional Information:

1. Pregnant women should not perform whitening. 2. Do not use tobacco products or eat during whitening. 3. White spots on the enamel that appear during whitening may disappear after a few days. 4. Stains existing on the teeth before whitening will not completely disappear. 5. Naturally, teeth are darker along the gum line and usually remain slightly darker. 6. A small percentage of patients experience sensitivity during whitening. If this happens, remove the trays and contact your dentist. 7. Regular dental check-ups and cleanings before and after whitening are important to maintain a healthy smile. 8. Some patients may notice temporary discomfort in the gums, lips, throat, or tongue. If any of these symptoms last for more than two days or gradually worsen, contact your dentist. These side effects typically disappear within 1-3 days after stopping the treatment. 9. Coffee, tobacco, and other products can cause tooth stains. If this happens, you can re-whiten your teeth using WHITE ANGEL. Foods and juices rich in citric acid may cause tooth sensitivity. 10. Old amalgam fillings may leave a gray color on your whitening trays. 11. Crowns, bridges, partial dentures, veneers, and composite fillings will not whiten. 12. Store the whitening kit away from sunlight and heat. Refrigeration is recommended, but do not freeze. 13. Keep the whitening gel out of reach of small children. 14. Do not swallow the gel. The product contains peroxide and may contain fluoride/potassium nitrate; swallowing large amounts can be harmful.

Note: This text is written in masculine form for convenience but is intended for both men and women.